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CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl

INDEX

SI No.	TOPICS	Pg No.
1.	What is Laboratory Manual	01
2.	Introduction	02-03
3.	Instruments used in Laboratory	4-14
4.	Maintenance of discipline cleanliness in the laboratory	15-16
5.	Exercise	17-26

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl


Ques:- What is Laboratory Manual ?

Ans:- Laboratory Manual is an exercise book containing various topics as allotted in the syllabus. This provides a prospective to students for practical work in the classrooms. By studying such manual a student becomes aware of the brief theoretical information, preparation of the sample material for study and observation, thereafter description and ultimately derivation of conclusion. The manuals are also provide with necessary diagrammes and sketches of the prepared sample as per observations made by the students.

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl

[Signature]
IQC
CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl

Introduction

2

→ The early botanical knowledge leads to the development and evolution of multidisciplinary subject as modern botany in recent years. The input from various sources of science and technology have tremendously enriched the subject Botany. Today this subject is a fusion of traditional and classical components with that of modern disciplines resulting in to a splendid natural science and technology have tremendously enriched the subject Botany. Today this subject is a fusion of traditional and classical components with that of modern disciplines resulting in to a splendid natural science. We have

CO-ORDINATOR

Vijaya Raju Govt. Girls PG College
Mogar-GWL

3

today various branches of Botany which have been enriched through inputs from extensive studies and researches, eg. Taxonomy / Systematics, Palaeobotany, Angiosperms Morphology, Ecology, Biodiversity and conservation, Palynology, Cytology and Cell Biology, Genetics, Plant Physiology and Biochemistry, Molecular Biology, Biotechnology, Microbiology, Phycology, Mycology and Plant Pathology, Lichenology, Bryology, Pteridology, Gymnosperms, Economic Botany, Pharmacognosy or Medical Botany, Ethnobotany etc.

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl


Instruments usually used in Botanical Lab. :-

1. Simple Microscope :

■ It consists of only one lens unit may even be an ordinary magnifying glass. Dissecting microscope is used either for dissecting the material or for less magnifications, i.e., only $5\times$, $10\times$ or rarely $20\times$. It is mainly used for taxonomic studies, embryo separation etc.

■ It consists of a basal foot and a limb. The 'stage', made up of a simple glass plate is attached to the limb. For the light adjustment purposes, a mirror is attached to the limb under the stage. Mirror can be moved vertically with the help of adjustment.

Compound Microscope

CO-ORDINATOR
 Vijaya Rajee Govt Girls PG College
 Morar-Gwl

At the tip of the limb is present a folded arm, on which a lens of definite magnification (5x, 10x etc.) is fitted. folded arm is moved to keep the lens in the desired position on the stage.

Light Microscope (compound) :

It is so named because it consists of two or more lens systems.

At the top is present the ocular lens. It can be turned around or may be removed. At the top of ocular lens is written 5x or 10x signifying the 5 times or 10 times magnification respectively.

Just below the ocular is a body tube, the bottom end of which contains a circular piece called objective lenses. Nose piece

6
The flat platform present below the objective is called **Stage**.

On the arm of the microscope are present two knobs called **Coars adjustment knob** and **fine adjustment knob**.

Out of the three objectives, the shortest is the low power objective. It has the largest lens but its magnifying power is least of the objective lenses. On the objective is high power of $10\times$ similar to ocular lens. It means if a $10\times$ ocular lens is used the magnification is $10 \times 10 = 100$ times.

The other objective is high power objective. Its magnification is equal to the number written on it multiplied by the power of ocular, i.e., $5\times$ or $10\times$ (objective $\times$ ocular).

⑦
The third objective is called oil immersion.

Generally, it contains a black band around the lower end. Use a drop of oil on the slide at the time of studying with the oil immersion objective. Its magnification can be estimated as $\text{ocular} \times \text{objective}$.

The use of oil (e.g., cedar-wood oil) is essential in order to keep the light rays properly aligned with the small objective.

Just below the stage is the condenser. Its function is to gather light from the mirror and to direct it to the objective lens. Condenser may be lowered or raised by a knob present on one side of the microscope beneath the stage.

Condenser contains a shutter called Iris diaphragm.

÷ Precautions for handling :-

1. Clean the ocular and objective lenses with lens paper / tissue paper, and do not remove them.
2. While studying an object, learn to keep one hand on the fine focus knob and focus continuously up and down.
3. While observing any kind of preparation, do not tilt the microscope.
4. Leave the low power objective in place after finishing all the observations.
5. To examine an object, always first use the low power and then the other objectives.
6. Never allow an objective lens to strike either the stage or a slide while focusing.
7. Use always the fine adjustment with high power objective.
8. Do not use oil immersion objective ~~without~~ oil.

CO-ORDINATOR

Vijaya Raje Govt Girls PG College
Morar-Gwl

Weighing Balance


Fig 3: Analytical Balance

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl

-:- Weighing Balance -:-

1

- Laboratory balance are defined by high levels of accuracy and precision in analytical testing and quantitative analysis.
- Used in a variety of applications, they are divided into five subcategories according to their readability:
 - Precision balance ($\geq 0.001\text{g}$)
 - Analytical balances (0.0001g)
 - Semi-micro balances (0.00001g)
 - Micro balances (0.000001g)
 - Ultra-micro balance (0.0000001g)

Precautions :-

- The instruments should be kept on a vibration free support.
- Every time use a clean spatula of appropriate size.
- Properly clean the surface after use.
- Always weigh sample after closing the weighing chamber doors.

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morang, Gwal

Single Glass Distillation Unit

Distilled water [DH_2O] is needed for different laboratory experiments such as culture media preparation, staining solutions preparation etc. Single glass distillation unit provides pyrogen free distilled water.

Precautions :-

Do not operate the machine water.

Always maintain the water level.

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl

-!- Autoclave -!-

- It is a pressure chamber used to carry out scientific process required elevated temperature and pressure in relation ambient.
- Used to sterilize equipment and culture media by subjecting them pressurized saturated steam at 121°C (250°F) for around 15-20 minutes.

Precautions :-

- Do not operate the machine without water.
- use only distilled water.
- Always maintain the water level above the heater.
- Never sealing containers; under pressure they pose an explosion risk.
- Never opening the door to the autoclave if there is at or below water running out the bottom clogged opening steam lines, equipment malfunction, or plugged drain may cause buildup for scalding water.

Volubility

- LCD Control Panel
- Ultraviolet lamp
- HEPA
- Primarily efficient filter
- Convenience receptacle
- Caster

Fig 8: Laminar Flow Cabinet

Laminar flow Cabinet

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl

-:- Laminar Flow Cabinet :-

12

A Laminar flow cabinet is a carefully enclosed bench designed to prevent contamination of semiconductor wafers, biological samples, or any particle sensitive material.

Air is drawn through a High-efficiency particulate air [HEPA] filter and blown in a very smooth, laminar flow towards the user.

Laminar flow cabinets have a UV-C germicidal lamp to sterilize the interior and contents before usage to prevent contamination of experiment.

Precautions :-

Germicidal lamp are usually kept on for 15 minutes to sterilize the interior.

No contact is to be made with a laminar flow hood during work.

Switch off UV lamp during use.

pH Meter

13

It is a scientific instrument that measures by the hydrogen-ion activity in water-based solution, indications, indication its acidity or alkalinity expressed as pH.

Precautions :-

pH electrode is fragile and should not be tempted to use it as a stirring glass rod when adjusting pH.

It needs to be calibrated daily with the help of standard buffer solutions.

Glass is rapidly attacked by hydrofluoric acid so never keep the electrode in HF solution as it lead to irreparable damage of electrode.

Adopt the suggested cleaning procedure after use
General cleaning $\rightarrow$ if the solution are fairly clear, soak in 0.1M HCl or 0.1M HNO_3 for about 30 min

• Bacter or deposit - soak in 1:10 dilution of laundry bleach for 10 min

• Protein deposit - soak in 1% pepsin in 0.1M HCl 5 min

oil or grease - soak in detergent or Methanol.

After every soaking the electrode should be rinsed with distilled water.

Incubator

14

This instrument used to grow and maintain microbiological cultures or cell cultures at desired temperature.

Precautions :-

- Place your unit properly.
- Clean the incubator regularly.
- Calibrate on a regular basis.

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl

Maintenance of discipline and cleanliness in Lab.

A list of do's and don'ts can be a harmful reminder safety issues. These are as follows-

Observe to :

15

Become familiar with experiments before coming to the laboratory.

Follow all laboratory instruction carefully.

Know the emergency and fire procedures.

Ask for clarification on any laboratory procedure that are not understood.

Wear safety glasses during practical where glassware, heat or chemicals are used.

Wear laboratory aprons.

Wear protective gloves when instructed to do so.

Keep all laboratory exits clear of debris.

Remove cultural spills and clean the area with 70% ethanol.

Entry in lab record book before using instruments.

Wash hands before leaving the laboratory.


Avoid to :-

- Begin laboratory work until instructed to do so.
- Do anything in the laboratory that is not understood.
- Engage in any horseplay in the laboratory.
- Perform unauthorized experiments.
- Work in isolation - a second person should always be nearby.
- Speak during work in Laminar Air flow / culture room.
- Touch your face when working with cultures.
- Light a flame source near flammable materials.
- Move a lit Bunsen burner.
- Taste or ingest any chemicals or plant material.
- Handle or consume food or drink in the laboratory.
- Pipet anything by Mouth.
- Pour chemicals back into dispensing containers.
- Leave a test tube stopper in place when heating the tubes.
- Drink from laboratory glassware.

[Gram-Positive Bacteria]

Fig 1.1: Gram staining of *Bacillus* sp. (Microscopic view).

[Gram-Negative Bacteria]

Gram Staining technique

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl

Gram staining technique

Requirements :

- 24 hours cultures of *Bacillus* sp. and *Escherichia coli*.
- Gram staining reagent :
 - Gram crystal violet solution.
 - Gram's iodine solution.
 - Safranin
- Ethanol (95%)
- Staining tray
- Wash bottle of distilled water
- Inoculating loop.
- Glass slides.
- Blotting paper.
- Bunsen burner / spirit lamp
- compound Microscope.

Procedure:

18

Make thin smears of one loop full inoculum of Bacterial samples on separate glass slide.

Let the smears air dry.

Heat fix air dry.

Cover each smear with crystal violet for 30 seconds.

Wash each slide with distilled water for a few seconds, using wash bottle.

Cover each smear with Gram's iodine solution for 60 seconds.

Wash off the iodine solution with distilled water.

Add Ethanol [95%] drop by drop, until no more colour flows from the smear. [The gram-positive bacteria are not affected while all gram-negative bacteria are completely decolorized.]

Wash the slides with distilled water and drain.

Apply safranin to smears for 30 seconds.

Wash with distilled water and blot dry with paper.

Let the stained air dry.


Observations :-

19

- Examine the slides microscopically using oil-immersion objective.
- Identify the gram reaction of both the culture and classify them.
- Make sketches for morphology of the cultures.
- Describe the morphology and arrangement of the cells.

Results :-


- Those bacteria that appear purple are referred to as gram-positive, those appearing pink are described gram-negative. In B.sp. the appear dark purple or blue in colour, thus it is a gram-positive bacterium whereas in E.coli the rods appear pink and is thus a gram-negative bacterium. To get reliable results, one should use cultures that are 18 to 24 hours old.

Blue-green Algae

Algae

[Signature]
CO-ORDINATOR
 Vijaya Raje Govt Girls PG College
 Morar-Gwl

From a
 Practical
Obj
 Preparation
 familiar
 synophy
Requ
 Compo
 Slide
 Needle
 10% G
 Cotton
 Vegetat
 Bunsen

EXERCISE : 2

20

Algae

From Algae one Nostoc is considered here for Practical work.

Objective :-

Preparation of Algae slide.

familiarity with its vegetative structure under cyanophyceae.

Requirements :-

Compound Microscope. • Blotting paper

Slide, Cover slip. • Wax

Needle and Brush

10% Glycerine


Cotton Blue solution [1%]

Vegetative part of Nostoc sp.

Vijaya Rajee Govt Girls PG College

Morar-Gw

Burner burner / spirit lamp

Method of Sample preparation :-

Place on the slide a little part of Nostoc by using 10% glycerin solution.

Spread the material with the help of needle and add a single drop of cotton blue [1%].


Wait for a while.

Take another slide, add few drops of lactophenol and place the spreaded material.

Place cover slip over the material carefully.

Care should be taken to absorb the excess amount of lactophenol without specimen blotting paper.

Place the prepared slide on turn and seal the margin of cover slip with melted with the help of a brush.

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl

Description and Identify characters :-

22

- Place the prepared slide under a compound microscope, draw the characteristics feature and observe the following features as identify character.

Thalli are present in the form of colony.

Ball like colony is enveloped by a gelatinous sheath.

Each colony contain thousand of straight or twisted filaments or Trichomes.

A trichomes is contorted and consists of many cells arranged in a beaded mannual.

Each cell is somewhat cylindrical or spherical shape.

There are present some large, spherical or cylindrical, colourless empty cell. present in the filaments, called heterocysts.

Some thick-walled cells of the filament become enlarged. These are called

÷ fungi ÷

from fungi one specimen, *Ascobolus* is considered here for practical work.

Objective :-

- Preparation of fungal slide.
- Familiarity with its fruit body structures containing Ascus and Ascospore under Ascomycete.

Requirements :-

- Compound Microscope.
- Bunsen burner/spirit lamp
- fruit body of *Ascobolus* sp.
- wax
- Slide, cover slip.
- Blade, Needle and Brush.
- Cotton blue solution [2%].
- Lactophenol [2%].
- Blotting paper.

EXERCISE : 3

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl

Method of Sample preparation

Make a thin vertical longitudinal section (VLS) of fruit body of the specimen with a sharp blade.

Transfer the section into the slide & add 1% cotton Blue solution.

wait for a while,

Place cover slip over the material carefully.

care should be taken to absorb the excess amount of stain by blotting paper.

Place the prepared slide on turn table & seal the margin of cover slip with melted wax with help of a brush.

Description and Identify characters

Place the prepared slide under a compound microscope, draw the characteristic features & observe the following features as identify characters.



Fruit body is Apothecium in nature.

Asci & Paraphyses are arranged alternately parallel.

Each ascus contain eight oval shaped ascospores.

Bryophyte

Marchantia

EXERCISE : 4

25

-:- Bryophyte -:-

From Bryophyte one specimen, *Marchantia* is considered here for practical work.

Objective :-

familiarity with gametophytic plant body. *Marchantia* Sq. under hepaticoids.

Requirements :-

Compound microscope

Permanent slide of antheridiophore of *Marchantia* Sq.

Description & Identify characters :-

Place the permanent slide under a compound microscope, draw the characteristic features & observe the following features as identify characters.

Antheridiophore consists of a long stalk and a lobed disc at apex.

The disc is usually eight lobed.

The disc consists of air chambers alternating with
hercoidal cavities, 26

Air chambers are more or less trigonular & open on upper surface by a pore called ostiole.

Antheridia in acropetal succession, the older near the centre & youngest at the margins.

CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morar-Gwl

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Under the Supervision of

Dr. PREETI GUPTA

(Assistant Professor, Botany Department)

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PROJECT WORK ON

"Study of Instruments used in Biotechnology Laboratory"

FOR THE PARTIAL FULFILMENT OF REQUIREMENT FOR

BACHELOR'S DEGREE

SESSION 2023-24

S.No.	NAME OF STUDENTS	CLASS	ROLL NUMBER
1 ^o	Kumkum Dantre	BSC 1 st year	Kumkum
2 ^o	Swati Sharma	BS.C 1 st yrs (Bio)	Sharma
3 ^o	Kajal Shrivastava	" " (Bio)	Kajal
4.	Priya Kain	" " (Biotech).	Priyakain
5.	Kajal Kushwah	" " (Biotech)	Kajal Kushwah
6.	Priyanshi Sharma	" " (Biotech)	Priyanshi
7.	Kavita Baghel	" (Botany)	Kavita
8.	Anjali Kaurav	" (Botany)	Anjali
9.	Kanchan	" " (Botany)	Kanchan
10.	Sangam Jayant	" " (Botany)	Sangam
11.	Poonam Kaurav	" " (Botany)	Poonam
12.	Arshi	" (Botany)	Arshi
13.	Kashish	(Ex) BSC 1 st yrs (Ex)	Kashish
14.	Ankita Tomar	FROM BS.C 1 st year	Ankita

DEPARTMENT OF BIOTECHNOLOGY

GOVERNMENT VIJAYARAJE POST GRADUATE COLLEGE

UNDER THE SUPERVISION OF

Mrs. Sarika Saxena (GUEST FACULTY LECTURER)

CO-ORDINATOR

Vijaya Raje Govt Girls PG College
Morar-Gwl


GOVERNMENT VIJAYARAJE GIRLS POST GRADUATE COLLEGE
JIWAJI UNIVERSITY

CERTIFICATE OF ORIGINALITY

This is to certify, that the Project Work Entitled "*Study of Instruments used in Biotechnology Laboratory*" submitted by us is an outcome of our independent and original work. We have duly acknowledged all the sources from which the ideas and extracts have been taken. The project is free from any plagiarism and has not been previously submitted elsewhere.

S.No.	NAME OF STUDENTS	CLASS	ROLL NUMBER	SIGNATURE
1.	Rumkum Dantre			Rumkum
2.	Swati Sharma			Swati
3.	Kajal Shrivastava			Kajal
4.	Kajal Kushwah			Kajal
5.	Priyanshi Sharma			Priyanshi
6.	Kavita Baghel			Kavita
7.	Anjali Kaurav			Anjali
8.	Kanchan			Kanchan
9.	Sangam			Sangam
10.	Poonam Kaurav			Poonam
11.	Arti			Arti
12.	Kashish			Kashish
13.	Priya Kain			Priya
14.	Ankita Tomar			Ankita

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APPROVAL LETTER

I hereby certify that the above mentioned reports are genuine reports of the Project Work done by the Students under my direction. This Project Work has been submitted to the college after my approval.

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supervisor

Vijayaraje Govt. Girls Post Graduate College

Morar, Gwalior (M.P.)

Department of Biotechnology
E-mail: vijayarajegirlscollege@gmail.com

Telephone No. 0751-2368329

CERTIFICATE OF COMPLETION

This is to certify that
.....listed.....students.....of.....B.Sc. I year.....(Botany
and Biotechnology).....
.....
.....
.....

have successfully completed their Project Work entitled *Study of Instruments used
in Biotechnology Laboratory* from Department of Biotechnology.

All the students are hard working, Deductive and results oriented. They all did
outstanding work during their tenure; we wish them a bright future.

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Enrollment No. $\Rightarrow$ JS23018909


Date $\Rightarrow$


CO-ORDINATOR
Vijaya Raje Govt Girls PG College
Morhar-Gw

Mo. No. $\Rightarrow$ 8827505805

INTRODUCTION

The microscope is a valuable instrument. There are many small objects or details of objects which cannot be seen by the unaided human eye. The microscope magnifies the image of such objects thus making them visible to human eye. Microscopes are used to ~~by~~ observe the shape of bacteria, fungi, parasites and host cells in various stained and unstained preparations. Microscope has the capacity to distinguish small gaps between two separate points which humans cannot distinguish. It is called its **resolving power or resolution**.

Microscope is used by scientists and healthcare professionals for many purposes, including diagnosis of infectious diseases, identification of **microorganisms** (microscopic organisms) in environmental samples (including food and water), and determination of the effect of pathogenic (disease-causing) microbes on human cells. This exercise will familiarize you with the microscopes.

A microscope (Greek: **mikron** = small and **skopos** = look)


HISTORICAL BACKGROUND

Topic _____ Date _____

- Hans and Zacharias Janssen -
- 1590, Dutch Eyeglass makers Inventors
- First compound microscopes (2 lenses), tube with lenses at each end.
- Robert Hooke -
- 1635-1703, English chemist, Mathematician Physicist and Inventor
- Compound microscope improvement
- Antonie van Leeuwenhoek -
- 1632-1723, Wine Assayer Surveyor, cloth Merchant, Mirror public official and Inventor
- Simple Microscope (1 lens)
- Father of Microscopy
- Ernst Ruska and Max Knoll -
- a physicist and an electrical engineer, from the University of Berlin.
- who created the first electron microscope in 1931.
- Manfred von Ardenne -
- was a German scientist and inventor
- who invented first scanning electron microscope.

VARIABLES USED IN MICROSCOPE

- Degree of enlargement

- ## 2. RESOLUTION -

- ### 3° NUMERICAL APERTURE (NA)

Designer Sheet 90

PARTS OF MICROSCOPE

1° Base :

It is a supporting stand, rests on the table and bears the weight of the microscope.

2° Arm (Limb) :

It is a curved, solid piece, movably adjusted with the base at the inclination joint. The arm holds the body tube and the stage.

3° Inclination joint :

It holds the body of the microscope with the base and permits inclination of the upper part of adjust to eye piece.

4° Body Tube :

It is a metal tube blackened inside. A revolving nose piece carrying the objectives is attached to the lower end. At the top of the body tube is screwed a draw tube which houses the eye piece.

5° Stage :

It is a sheet of squarish or rectangular platform with 10 to 12 cm sides and having a circular hole at the centre. It is made up of highly polished metal or bakelite and provided with two spring clips for holding the slide or sliding tray. The stage is firmly secured with the arm.

6° Coarse adjustment -

It is attached to the arm. The adjustment has a large head. Clockwise turn of the pinion head moves the body tube downwards and anticlockwise turn moves the body tube upward.

7° Fine adjustment -

It is also attached to the arm, parallel to the coarse adjustment. It operates in the way as in the coarse adjustment, but the movement of the body tube is only slight with the turning of the pinion. Usually ten turns of the fine adjustment pinion are equal to one turn of that of the coarse adjustment.

8° Mirror

A planconcave round mirror is fitted below the stage at some distance from it. In artificial light the plain surface is used. In skylight the concave surface is used, as more light is reflected to the object from this surface. It required the mirror can be taken out from the hole in the base of the microscope and a sub-stage lamp fitted.

9° Illuminating apparatus (Condenser) -

It is fitted below the stage and provided with an iris. The opening of the diaphragm can be reduced or even closed completely. The arrangement is to regulate the amount of light passing through the object.

10. Diaphragm or Iris :-





The diaphragm or Iris is located under the stage and is an apparatus that can be adjusted to vary the intensity and size, of the cone of light that is projected through the slide. As there is no set rule on which setting depends on the transparency of the specimen and the degree of contrast you desire in your image.

11. Rack and Pinion -

To bring the object under focus, the rack and pinion are either attached to the body tube or the stage.


12. Draw Tube :-

The upper end of the body tube has a small fixed tube which is known as the draw tube. The main function of the draw tube is to hold the ocular lens.

THE EFFECTS OF IMMERSION OIL ON RESOLUTION

[Signature]

CO-ORDINATOR

Vijaya Rajee Govt Girls PG College
Morar-Gwl

TYPES OF MICROSCOPE

Topic _____ Date _____

• Light microscope : Use ~~light~~ sunlight or artificial light.

(A) Bright field microscope

(B) Dark field microscope

(C) Phase contrast microscope

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• Electron microscope :

(1) Transmission electron microscope

(2) Scanning electron microscope

A° Bright field microscope -

→ Produces a dark image against a lighter background.

→ Has several objective lenses.

→ 2 types : Simple
Compound

• Simple -

Contain a single magnifying lens.


• Compound -

Series of lenses for magnification

Light passes through specimen into objective lens

oil immersion lens increases resolution.

IMAGE FORMED BY BRIGHT FIELD MICROSCOPY

Light (optical) microscope

using visible light to generate images

Simple microscope

Leeuwenhoek microscope (1600s)

having more than one lens

Compound microscope

Stereo microscope

Digital microscope

Monocular microscope

Binocular microscope

stereoscopic view

Dissecting microscope

digital imaging

USB digital microscope

viewing fluorescent objects

viewing specimens from the bottom

[Signature]

CO-ORDINATOR

Vijaya Raje Govt Girls PG College
Morar-Gwl

Have one or two ocular lenses.

Resolution = 200nm .

Advantages -

- Used to view live or stained cells.
- Simple setup with very little preparation required.

Disadvantages -

- Biological specimens are often of low contrast and need to be stained.
- Staining may destroy or introduce artifacts.
- Resolution is limited to 200nm .

B. Dark field microscope -

- Produces a bright image of the object against a dark background.
- Optical system to enhance the contrast of unstained bodies.
- Specimen appears gleaming bright against dark background.

Requisites For Dark field microscopy -

- Dark ground condenser
- High intensity lamp
- Funnel stop

CARE
OF
MICROSCOPE

Topic _____

Date _____

- The following rules, cautions and maintenance hints will help keep your microscope in good operating condition -

1° Use both hands when carrying the microscope one firmly grasping the arm of the microscope, the other beneath the base. Avoid jerking your microscope.

2° Never touch the lenses. If the lenses become dirty wipe them gently with lens tissue.

3° If blurred specks appear in the field of view this may be due to lint or smears on the eyepiece. If the specks move while rotating the eyepiece the dust is on the eyepiece and cleaning the outer lens of the eyepiece is in order. If the quality of image is improved by changing objective lenses, clean the objective lens with lens paper.

4° Never leave a slide on the microscope when it is not in use.

5° Always remove oil from the oil-immersion objective lens after its use. If by accident

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
CO-ORDINATOR

Vijaya Raj Govt Girls PG College
Morar-Gwl

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1. introduction.
2. principle
3. construction
4. working.

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
B.S.C 1st YEAR


DEFINITION

A water Bath is a piece of Scientific Equipment used to maintain a steady temperature for a prolonged time when incubating samples. ⁹⁹

A water bath is preferred over an open flame when heating flammable substances. It employed to enable some chemical rxn at high temperature

Principle of WATER BATH.....

The sensor converts the temperature of the water into a resistance value, which is then amplified. This produces the control signal, which effectively regulates the average heating power of the electric heating tube and keeps the water at a constant temperature.

DFY structure

Big container filled with warm water makes
 a hot water bath. It has different design
 combination, sizes, proportion. A laboratory water
 bath's container capacity ranges from 12 to 32
 litres for a standard model & 50 to 100 litres
 for a large water bath. Most water baths feature
 a digital or analog interface that users may use
 to set the desired temperature. In contrast, some water
 baths include a reader that uses a current to regulate

he temperature. A shaker is another feature of water bath that is useful for homogenizing the test solution.

It is mostly utilized in clinical & microbiology labs, university labs, Environmental research even food technology for warming reagents, Sample mixing, corrosion tests, bacteriological investigations, among other applications. A water bath is tool that needs to additional supporting instruments to functions.

An incubator, a micropipette, a centrifuge, a laboratory are typically used in tests involving a water bath.

Certificate

Name Purva Jain.

Class : B.Sc 1st Year.

Roll No :

Exam No. :

School _____

*This is certified to be the bonafide work of the student in the
..... Laboratory during the academic
year 20 / 20 .*

*No. of practicals certified out of in the
subject of*

.....
Teacher In-charge

.....
Examiner's Signature

.....
Principal

Date:

School Rubber Stamp

(N.B. The candidate is expected to retain his/her journal till he/she passes in the subject)

Date _____

Expt. No. _____

Page No. _____

Field Project

Incubator

Teacher's Signature _____

Incubator

Definition.

An incubator in microbiology is an insulated and enclosed device that provides an optimal condition of temperature humidity and other environmental conditions required for the growth of organism.

An incubator is a piece of vital laboratory equipment necessary for cultivating microorganisms under artificial conditions.

An incubator can be used to cultivate both unicellular and multicellular organisms.

A microbial incubator is made up of various units some of which are :-

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